

Four Easy-to-Spot Differences Between EXPAREL and DIPRIVAN® (Propofol)

Vial size and shape

Short and wide

1

DIPRIVAN

Cap size, color, mechanism

Wide cap; yellow and blue

Metal ring and stopper remain intact when cap is flipped off

2



Label colors, font, directives

Black letters on white background

Bold black lettering: For intravenous administration

Bold red lettering: Shake well before using

3

Dosage

1%; 200 mg/20 mL (10 mg/mL)

4

EXPAREL

1

Vial size and shape

Tall and narrow; shortened neck space

2



3

Label colors, font, directives

Teal letters on white background

Teal bar with bold white lettering: For infiltration and interscalene brachial plexus nerve block. Not for any other route of administration

Gray banner surrounding NDC and 1.3% concentration

4

Dosage

1.3%; 266 mg/20 mL (13.3 mg/mL)

Please see full Prescribing Information for EXPAREL, available at www.EXPAREL.com.
For more information, please visit www.EXPAREL.com or call 1-855-793-9727.